

DS PrimeTaper implant system: reliable, efficient, seamless

DS PrimeTaper implants enable optimized implant treatments, delivering desired primary stability at placement, which enables immediate restorations. Bone care concept promotes reliable function and esthetics, with long-term marginal bone maintenance. This new system provides an efficient handling experience and is supported by a seamless, fully integrated digital implant workflow.

	Description
Implant	Tapered design with variable thread shapes OsseoSpeed surface Soft Tissue Chamber EV Connection with conical seal design
Sizes	Diameter from 3.6 to 5.4 mm Length from 6.5 to 17 mm
Bone Care Concept	Progressive design for initial stability Linear torque build-up throughout the implant installation Stable and fast implant installation procedures thanks to double thread design Pronounced flute design with high cutting capacity MicroThread for cortical bone maintenance
Primary and long-term stability	Progressive tapered design provides high primary stability allowing for immediate loading and temporization. Secondary stability and long-term bone preservation, ensured through proven features derived from Astra Tech Implant System EV
Efficient handling	Faster implant installation due to the double thread design Easy drilling protocol and logical tray layout Predictable abutment installation due to one-position-only connection New Implant Driver EV for secure and consistent pick-up and insertion of the implant

Press Contact

Guenther Schmidhuber
 Director Corporate Communications
 Sirona Straße 1
 5071 Wals bei Salzburg, Austria
 T +43 (0) 662 2450-566
 F +43 (0) 662 2450-540
guenther.schmidhuber@dentsplysirona.com

Allegra Cau
 Daniel J. Edelman, Holdings Inc.
 250 Hudson Street, New York
 NY 10013, USA
 T ++1 212-277-3746
dentsplysirona@edelman.com
www.edelman.com

Marion Par-Weixlberger
 Vice President Public Relations &
 Corporate Communications
 Sirona Straße 1
 5071 Wals bei Salzburg, Austria
 T +43 (0) 662 2450-588
 F +43 (0) 662 2450-540
marion.par-weixlberger@dentsplysirona.com

About Dentsply Sirona:

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with over a century of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better and safer dental care. Dentsply Sirona's headquarter is located in Charlotte, North Carolina. The company's shares are listed in the United States on NASDAQ under the symbol XRAY.

Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

Longevity and natural esthetics	<p>Comprehensive robust prosthetic solutions enable healthy soft tissue for long-lasting esthetics</p> <p>Efficient restorative workflows for all indications and treatment protocols</p> <p>Tight and stable fit of all restorative components, minimizing micromovements and microleakage</p>
---------------------------------	---

