

The Intego Treatment Center: optimally designed for preventive care and ready for the future

As part of its commitment to better, safer, faster dentistry, Dentsply Sirona is pleased to present another unique integrated solution. The Cavitron Ultrasonic Scaling System with Steri-Mate 360 is now fully integrated in the Intego Treatment Center.

York, August 1st 2017 – General scaling and periodontal debridement applications are an integral service offered by dental practices. In order to help practitioners to best set up their treatment rooms, Dentsply Sirona has expanded its existing Treatment Centers product range with an integrated preventive care concept. The Cavitron built-in ultrasonic scaling system is now available with the Intego line of Treatment Centers for fast and effective removal of even the most tenacious deposits. The Intego and Intego Pro Ambidextrous models are available in the US and offer modern dental practice everything necessary for professional dental hygiene, dental aesthetics and dental care.

Flexibility through individually configurable equipment

Intego supports a variety of treatment concepts through two different variants: Intego TS is a traditional, over the patient delivery with hanging hoses. The dentist element has a wide range of movement and the adjustable tray holder makes for a convenient work environment. Intego CS is a continental delivery system with whip arms and offers the widest range of treatment positions around the oral cavity, whether working alone or as a team.

Intego is additionally configurable with special packages and practical equipment options. For optimal patient communication, the unit can be equipped with a 22-inch full-HD monitor, and an autofocus intraoral camera.

Integrated design ideal for hygiene procedures

Intego is an excellent choice for hygiene procedures. Configuration is ergonomically and functionally adapted to preventive care treatments thanks to features such as a double articulating headrest, external device connection, multifunctional syringe, maintenance-free electric motors, LED operating lights and, coming soon, the integration of Cavitron. The Intego Pro Ambidextrous model features an integrated infection prevention system and automated cleaning functions. Ergomotion provides a synchronized movement of the patient's spine with no readjustment of the headrest. Changing from right-handed to left-handed operation can be easily done in 15 seconds between the treatments of patients. The Carl or Hugo stools are available as ergonomic add-ons.

Press Contact

Marion Par-Weixlberger
Director Corporate Communications
and Public Relations
Sirona Straße 1
5071 Wals bei Salzburg, Austria
T +43 (0) 662 2450-588
F +43 (0) 662 2450-540
[marion.par-
weixlberger@dentsplysirona.com](mailto:marion.par-weixlberger@dentsplysirona.com)

Kerstin Schicha
Edelman.ergo
Charlottenstraße 68
10117 Berlin
T +49 (0)30 201805-70
kerstin.schicha@edelmanergo.com
www.edelmanergo.com

ABOUT DENTSPLY SIRONA:

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 130-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company, Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY. Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

Ready for the future thanks to innovative design and digital interfaces

Developed at the Dentsply Sirona Innovation Center in Bensheim, Germany, Intego Treatment Centers were created as high-quality and reliable units, designed for the technical requirements of the future. A network interface in the patient chair, and USB connections in the doctor's element enable digital updates as well as the use of additional applications. New system solutions can be accessed and controlled via the comprehensive user interface on the EasyTouch touchscreen.

The Intego Treatment Centers combine modern, award-winning design available in many color variations, with compact dimensions, and a holistic infection prevention concept. Intego is optimally equipped for every clinical need while also being optimally prepared for long-term use into the future.

ILLUSTRATIONS

are available for [> download](#) on the website.



*Fig 1: The Caviron Ultrasonic System with Steri-Mate 360
- integrated with the Intego line of Treatment Centers.*

